



Stop Hém[®]

FR Pansement

HEMOSTASE ET CICATRISATION

12 PANSEMENTS :
6X (72 X 25 mm)
6X (72 X 34 mm)

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12 BANDAGES :
6X (72 X 25 mm)
6X (72 X 34 mm)

1. COMPOSITION / PROPERTIES

The sterile Stop Hemo® Bandage is a bioactive pad of calcium alginate fibres (from marine algae) on an adhesive support with high skin tolerance.

- The pad releases its calcium to accelerate platelet aggregation, stop bleeding quickly and promote healing.
- The calcium alginate fibres form a gel in contact with the blood, also ensures painless removal.

2. INDICATIONS

To stop the bleeding and promote the healing of superficial wounds: cuts, open blisters, small wounds. Stop Hemo® Bandage may be used in children, adults and pregnant women, and is recommended for patients with congenital (Rendu-Osler, von Willebrand, haemophilia) or acquired (under antiplatelet agents, anti-vitamin K or heparin) coagulation disorders. Persistent bleeding requires a medical consultation.

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Made in France by
Les Laboratoires Brothier

3. USAGE

Clean the wound with soap and water or physiological serum (NaCl 0,9%). Apply Stop Hemo® Bandage. *Fibers potentially remaining after Stop Hemo's removal shall be eliminated naturally.*

4. CONTRAINDICATIONS AND SIDE EFFECTS

No contraindication known to date. Stop Hemo® Bandage contains trace amounts of surfactant and its use may cause hypersensitivity reactions in some sensitised subjects.

5. INTERACTIONS

Stop Hemo® Bandage must not be used with alkaline solutions, e.g. sodium hypochlorite or solutions containing chlorhexidine, benzalkonium chloride or copper salts. If such products are used on the wound, thoroughly rinse with physiological serum (NaCl 0,9%) before applying Stop Hemo®.

6. PRECAUTIONS FOR USE

Do not reuse, because of the risk of microbial contamination and loss of performance. Always refer to the expiry date.

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